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[1] [http://www.thelancet.com/journals/lancet/article/PIIS0140-6736\(06\)69267-2/fulltext](http://www.thelancet.com/journals/lancet/article/PIIS0140-6736(06)69267-2/fulltext)

[2] Siew, A. "Pulmonary Drug Delivery – Particle Engineering for Inhaled Therapeutics", *Pharmaceutical Technology*, February 2014, pp.34-38.

[3] <http://conf.montreal-protocol.org/meeting/oewg/oewg-32/ngo-publications/NGO%20Publications/IPAC%20Position%20Statement%2032%20OEWG.pdf>

[4] Purewal TS. Test methods for inhalers to check performance under use and un-intentional mis-use conditions. *Pharmeuropa*. 2002;14(3).

[5] CPMP (2006) "Guideline on the Pharmaceutical Quality of Inhalation and Nasal Products" to be read in conjunction with –

CPMP (2004) "Points to consider on the Requirements for Clinical Documentation for Orally Inhaled Products (OIP)"

[6] Ph. Eur. Chapter 2.9.44 Preparations for nebulisation: Characterisation.

[7] USP 1601 Products for nebulization: Characterization.

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[9] M. Copley and P. Kippax. "From actuation to deposition: particle sizing techniques for characterizing nasal drug delivery systems", *Inhalation*, April 2012.

[10] J.P.Mitchell and M.W. Nagel (2004). Particle size analysis of aerosols from medicinal inhalers. *KONA Powder and Particle*, No. 22.

[11] USP 31 Chapter <601> Aerosols, Nasal Sprays, Metered Dose Inhalers and Dry Powder Inhalers; USP 31 Chapter <905> Uniformity of Dosage Units; USP 31 Chapter <1151> Pharmaceutical Dosage Forms (Aerosols – Inhalations)

[12] Ph. Eur. 6.0 Chapter 2.9.18 Preparations for Inhalation: Aerodynamic Assessment of Fine Particles; Ph. Eur. 6.0 Monograph 0671 Preparations for Inhalation; Ph. Eur. 6.0 Monograph 0523 Pressurised Pharmaceutical Preparations; Ph. Eur. 6.0 Monograph 0676 Nasal Preparations

[13] CDER (1998). "Metered Dose Inhaler (MDI) & Dry Powder Inhaler (DPI) Drug Products", *Chemistry, Manufacturing and Controls Documentation – Draft*

[14] M. Copley. "Cascade impactors – Theory, design and practical information for optimal

testing", *Inhalation*, February 2008, Vol 2, No 1, pp. 19 - 23

[15] Mitchell, J. P. et al (2003). Cascade Impactors for the Size Characterization of Aerosols from Medicinal Inhalers: Their Uses and Limitations. *Journal of Aerosol Medicine*, Volume 16, Number 4, pp. 341-377.

[16] M. Copley. "Cascade impactors enter a new decade", *Manufacturing Chemist*, May 2007, pp. 45-47