

Label Printer

As a space-efficient alternative to traditional A4 desktop printers, the Label Printer provides easy set-up, intuitive operation, and high-quality report printing, improving laboratory productivity. Printing all critical test parameters onto an adhesive report label, the Label Printer helps streamline laboratory data processing. With a simple interface and compact benchtop design, the Label Printer seamlessly integrates with compatible Copley products for efficient data reporting.

Key Features

-  Suitable for high-volume printing
-  Minimal maintenance required
-  High resolution adhesive report labels
-  Simple user interface
-  Compatible with Direct Thermal or Thermal Transfer Printer Labels*
-  Ideal for space-constrained labs




Label Printer shown here with the Tapped Density Tester JV 210i

Label Printer

Cat. No.	Description		
9765	Label Printer	9719	Thermal Transfer Printer Labels
9716	Direct Thermal Printer Labels	9725	Thermal Transfer Printer Ribbon

* The Label Printer (9765) is supplied with Direct Thermal Labels as standard. Fade resistant Thermal Transfer Printer Labels (9719) are also available for increased durability, (Thermal Transfer Printer Ribbon (9725) also required).

Inhaler Testing



Vertus® III, Vertus III+



DecaVertus® III



Critical Flow Controller
TPK 100i, TPK 100i-R
(firmware version 1.3.0 and above)



Sample Preparation Unit
SPU 200i
(firmware version 1.3.0 and above)




Breathing Simulator
BRS 100i
(firmware version 1.3.0 and above)



Breath Actuation Controller
BAC 100i, BAC 100i-R
(firmware version 1.3.0 and above)



Gentle Rocker™
GR 200i
(firmware version 1.3.0 and above)



DUSA Shaker™
DTS 100i
(firmware version 1.1.0 and above)

Pharmaceutical Testing




Tapped Density
JV 110i, JV 210i



Tablet Dissolution Tester
DIS 600i, DIS 800i
(firmware version 1.2.0 and above)



Tablet Hardness Tester
TBF 100i
(firmware version 1.1.0 and above)



Tablet Disintegration Tester
DTG 100i, DTG 200i, DTG 300i,
DTG 400i, DTG 200i-IS
(firmware version 1.2.0 and above)



Tablet Friability Tester
FRV 110i, FRV 210i